



Lucitone Digital Fit™ Denture Disc

The Lucitone Digital Fit™ Denture Disc enables the dental lab to produce a Dentsply Sirona Digital Denture in house. The lab can design an upper or lower digital denture arch, mill it, bond it with premium IPN 3D™ Digital Denture Teeth, polish and go.

Lucitone Digital Fit Denture Base Discs:

- Produce digital dentures recognized for precision & better fit.
- Maintain dimensional stability under high-speed milling.
- Offer a cleaner, safer way to fabricate dentures in your lab.
- 30 hours of bur life.

Lucitone is the industry benchmark for high strength, high-impact resistant denture base acrylics with over 40 years of clinical success.

- Premium choice for both tissue and implant-supported full and partial dentures offering high-impact resistance and superb flexural strength.
- Lifelike aesthetic and balanced translucency replicate natural gingiva.
- Available in seven shades with and without veins.

Dentsply Sirona Digital Dentures can be designed using inLab SW 20.0 or other design software.

- For more information, contact your Dentsply Sirona Lab Sales Representative.



Lucitone Digital Fit Denture Discs 98.5mm diameter for use in 5-axis milling machines

REF#	Size	Shade	Characterized	Application
D906110	20mm	Original	Veined	Denture Base Plate
D906114		Light	Veined	
D906118		Light Reddish Pink	Veined	
D906122		Dark Pink	Veined	
D906126		Original Opaque	Veined	
D906130		Pink Opaque	Non-Veined	
D906134		Pink	Non-Veined	
D906111	25mm	Original	Veined	Full Denture Arch Upper/Lower*
D906115		Light	Veined	
D906119		Light Reddish Pink	Veined	
D906123		Dark Pink	Veined	
D906127		Original Opaque	Veined	
D906131		Pink Opaque	Non-Veined	
D906135		Pink	Non-Veined	
D906112	30mm	Original	Veined	Full Denture Arch Upper/Lower*
D906116		Light	Veined	
D906120		Light Reddish Pink	Veined	
D906124		Dark Pink	Veined	
D906128		Original Opaque	Veined	
D906132		Pink Opaque	Non-Veined	
D906136		Pink	Non-Veined	
D906113	35mm	Original	Veined	Full Denture Arch Upper/Lower*
D906117		Light	Veined	
D906121		Light Reddish Pink	Veined	
D906125		Dark Pink	Veined	
D906129		Original Opaque	Veined	
D906133		Pink Opaque	Non-Veined	
D906137		Pink	Non-Veined	

* Select size based on denture design height

Digital Denture Workflow



REF#	Description	Application
D33135	Lucitone Family Shade Guide	Shade Selection

REF#	Description	Application
D6686286	Adapter for inLab MC X5, 3/PK	Disc Adapters

Dentsply Sirona | dentsplysirona.com

© 2020 Dentsply Sirona Inc. All rights reserved. Form No. DP-0000530 Rev. 1 (09/2020)