

Zirkon gesintert Zirconia sintered	Zirkon ungesintert Zirconia un sintered	Keramik Ceramic	Komposit Composite	Prothesen- kunststoff / Acryl Denture acrylic	NEW-Legierungen Non-precious alloys	Edelmetall / Amalgam Precious metals / Amalgam	Titan Titanium
	■						

- ☑ Speziell zur Entfernung der Frässtege
- ☑ Stark reduziertes Ausschußrisiko
- ☑ Ergonomie mit Zahntechnikern abgestimmt
- ☑ Specially designed for removal of milling bars
- ☑ Significantly reduced scrap risk
- ☑ Ergonomic developed with dental technicians



labor

BWG | M/F Artikel/Item

Körnung mittel/fein
Grit medium/fine

BWG164104MF



Ø 10 x 22,7

104-HP

8.000 / 20.000

Maße | Sizes [mm]

Schaft | Shank

○ optimal / maximum [min.-1]

ISO ISO-Nummer Seite 46 – 48
ISO-Nummer page 46 – 48

Der einzigartige Weisslingspolierer ist in enger Zusammenarbeit mit Fachleuten und Laboren direkt in der praktischen Anwendung entstanden.

Die Besonderheiten sind die auf die Aufgabe abgestimmte spezielle Form des Polierers, sowie die zwei assemblierten Schleifmittelzonen. Die graue Zone reduziert schnell und kontrolliert die Frässtege. Mit der farbpigmentlosen weissen Zone wird der präzise Feinschliff gearbeitet, ganz ohne die Gefahr von Oberflächeneinfärbungen und Beschädigungen des teuren Zahnkörpers.

Vorteil: Präzises, effizientes und angenehmes Arbeiten bei stark reduzierter Beschädigungsgefahr.

The unique Weissling polisher was developed in close collaboration with experts and laboratories working directly with such practical applications.

The special features include the precise form of the polisher, which is perfectly tailored to the task and its two assembled abrasive material zones. The grey zone reduces quickly and controlled the milling bars. The pigment-free white zone is used for precise fine-polishing, with no risk of surface colouring or damage to the precious working part.

The result: precise, efficient and pleasant work with a significantly reduced risk of damage.



Zone Grau: schnelle Reduktion der Frässtege.

Grey zone: quick reduction of the milling bars.



Weißer Zone: Schritt zwei – derselbe Polierer. Risikofreie Konturierung durch geminderte Abtragsleistung.

White Zone: Step two – same polisher. Riskless contouring through reduced removal.



Zwei in Einem: schnelle, präzise und sichere Bearbeitung.

Two in One: fast, accurate and reliable processing.